


**Quantitative Cellular and Molecular Biology
Laboratory
Computational Biology Department
Comp Bio 02-261**

PCR Primer Design Lab Introduction

Gel Electrophoresis

- Experimental method to determine distribution of DNA strand sizes in DNA sample.
- More details in the next lecture...





DNA Notation (primary and secondary structure)

ssDNA = 5' -ACTGCGATAGACGATGTCCGGATGACA-3'

← Shows sequence


dsDNA = 5' -ACTGCGATAGACGATGTCCGGATGACA-3'
3' -TGACGCTATCTGCTACAGGCCTACTGT-5'

← Shows sequence
and pairing



dsDNA =
5' -  -3'
3' -  -5'

← Shows pairing

Polymerase Chain Reaction



Polymerase Chain Reaction



Millions of copies!



2nd Cycle

Video URL

<https://www.youtube.com/watch?v=YJKYSIJREIc>

Tasks for Computational Lab

1. Generate features to allow prediction of primer melting points
2. Implement function for predicting PCR products
3. Design primers for PCR reaction to identify three types of DNA

Task 1 – Primer Melting Point Prediction

Features: Numerical descriptors of an object

Design ***features*** to help predict the melting point for a primer.
Implement your feature calculation methods.

Assess with N-fold cross-validation using a RandomForest regressor model for generating predictions.


Regressor/regression – predict continuous value

Classifier – predict discrete class

How to design features for predicting melting point?

http://www.premierbiosoft.com/tech_notes/PCR_Primer_Design.html

Assessing Accuracy of Predictions



$$R^2 = 1 - \frac{SS_{res}}{SS_{tot}}$$

R^2 values closer to 1.0 are better.

Task 2 – Predict PCR Products

- Given a sequence and a pair of primers, predict whether or not there will be a product. If so, predict the resulting product sequence (upper strand).
- Important Primer Pair Characteristics:
 - Reverse Primer
 - Binds to upper strand
 - Reverse compliment of binding site on upper strand
 - $T_m \sim 60^{\circ}\text{C} \pm 2.0$
 - 18-35 bases long
 - Forward Primer
 - Binds to lower strand within 1000 bases upstream of the reverse primer binding location
 - Reverse compliment of binding site on lower strand
 - $T_m \sim 60^{\circ}\text{C} \pm 2.0$
 - 18-35 bases long

Task 2 – Predict PCR Products

- Given a sequence and a pair of primers, predict whether or not there will be a product. If so, predict the resulting product.
- Important Primer Pair Characteristics:
 - Reverse Primer
 - Binds to upper strand
 - Reverse complement of binding site on upper strand
 - $T_m \sim 60^{\circ}\text{C} \pm 2.0$
 - 18-25 bases long
 - Forward Primer
 - Binds to lower strand within 1000 bases upstream of the reverse
 - Reverse complement of binding site on lower strand
 - $T_m \sim 60^{\circ}\text{C} \pm 2.0$
 - 18-25 bases long



Use Task 1!

Task 2 – Predict PCR Products

- Given a sequence and a pair of primers, predict whether or not there will be a product. If so, predict the resulting product.
- Important Primer Pair Characteristics:
 - Reverse Primer
 - Binds to upper strand
 - **Reverse compliment of binding site on upper strand**
 - $T_m \sim 60^{\circ}\text{C} \pm 2.0$
 - 18-25 bases long
 - Forward Primer
 - Binds to lower strand within 1000 bases upstream of the reverse primer
 - **Reverse compliment of binding site on lower strand**
 - $T_m \sim 60^{\circ}\text{C} \pm 2.0$
 - 18-25 bases long

Reverse Complement:
ACTG -> CAGT

Complementary Base Pairs:
A <-> T
G <-> C

Task 2 – Predict PCR Products

- Given a sequence and a pair of primers, predict whether or not there will be a product. If so, predict the resulting product.
- Important Primer Pair Characteristics:
 - Reverse Primer
 - **Binds to upper strand**
 - Reverse compliment of binding site on upper strand
 - $T_m \sim 60^{\circ}\text{C} \pm 2.0$
 - 18-35 bases long
 - Forward Primer
 - **Binds to lower strand within 1000 bases upstream of the reverse primer binding location**
 - Reverse compliment of binding site on lower strand
 - $T_m \sim 60^{\circ}\text{C} \pm 2.0$
 - 18-25 bases long

How do we determine binding?

Task 2 – Predict PCR Products

- Given a sequence and a pair of primers, predict whether or not there will be a product. If so, predict the resulting product.
- Important Primer Pair Characteristics:
 - Reverse Primer
 - **Binds to upper strand**
 - Reverse complement of binding site on upper strand
 - $T_m \sim 60^{\circ}\text{C} \pm 1.5$
 - 18-25 bases long
 - Forward Primer
 - **Binds to lower strand within 1000 bases upstream of the reverse primer binding location**
 - Reverse complement of binding site on lower strand
 - $T_m \sim 60^{\circ}\text{C} \pm 1.5$
 - 18-25 bases long

How do we determine
binding?

(Local) Sequence Alignment!

Task 2 – Predict PCR Products (Alignment)

```
> alignment.local_align("ACTG", "ACTG", print_output = True)
```

Scoring: match = 10; mismatch = -5; gap_start = 0; gap_extend = -7

A matrix =

	*	A	C	T	G
*	0	0	0	0	0
A	0	10	3	0	0
C	0	3	20	13	6
T	0	0	13	30	23
G	0	0	6	23	40

Optimal Score = 40

Max location in matrix = (4, 4)

Best Alignment:
ACTG
ACTG

Task 2 – Predict PCR Products (Alignment)

```
> alignment.local_align("ACTGACTGACTG", "ACTG", print_output = True)
```

Scoring: match = 10; mismatch = -5; gap_start = 0; gap_extend = -7

A matrix =

	*	A	C	T	G	A	C	T	G	A	C	T	G
*	0	0	0	0	0	0	0	0	0	0	0	0	0
A	0	10	3	0	0	10	3	0	0	10	3	0	0
C	0	3	20	13	6	3	20	13	6	3	20	13	6
T	0	0	13	30	23	16	13	30	23	16	13	30	23
G	0	0	6	23	40	33	26	23	40	33	26	23	40

Optimal Score = 40

Max location in matrix = (12, 4)

Multiple
Best
Alignments

Task 2 – Predict PCR Products (Alignment)

```
> alignment.local_align("AGTCACTGGCTT", "ACTG", print_output = True)
```

Scoring: match = 10; mismatch = -5; gap_start = 0; gap_extend = -7

A matrix =

	*	A	G	T	C	A	C	T	G	G	C	T	T
*	0	0	0	0	0	0	0	0	0	0	0	0	0
A	0	10	3	0	0	10	3	0	0	0	0	0	0
C	0	3	5	0	10	3	20	13	6	0	10	3	0
T	0	0	0	15	8	5	13	30	23	16	9	20	13
G	0	0	10	8	10	3	6	23	40	33	26	19	15

Optimal Score = 40

Max location in matrix = (8, 4)

Position in String 1
of the last character
in optimal
alignment

Position in String 2
of the last character
in optimal
alignment

Best Alignment:
----ACTG----
ACTG

Best Score:
40/40
Best score possible
for alignment of 4
characters.

Binding defined by 80%
+ alignment.

Local Alignment Function

```
def local_align(x, y, score=ScoreParam(10, -5, -7), print_output = False):
```

x = sequence 1

y = sequence 2

score = Score Parameter (match = +10, mismatch = -5, gap = -7)
(optional)

print_output = binary indicating whether or not you want pretty
output printed from alignment
(optional)

Task 2 – Predict PCR Products

- Given a sequence and a pair of primers, predict whether or not there will be a product. If so, predict the resulting product.
- Important Primer Pair Characteristics:
 - Reverse Primer
 - **Binds to upper strand (80%+ alignment)**
 - Reverse complement of binding site on upper strand
 - $T_m \sim 60^{\circ}\text{C} \pm 2.0$
 - 18-30 bases long
 - Forward Primer
 - **Binds to lower strand within 1000 bases upstream of the reverse primer binding location (80%+ alignment)**
 - Reverse complement of binding site on lower strand
 - $T_m \sim 60^{\circ}\text{C} \pm 2.0$
 - 18-30 bases long

How do we determine
binding?



(Local) Sequence Alignment!

Polymerase Chain Reaction Quiz!

5'	ACTGACGATACGATAGGCTACGAGCTCAGCGACTCATACG	3'	DNA Sequence
		3' GCTGAGTA 5'	Reverse Primer
5'	GACGATAC	3'	Forward Primer
3'	TGACTGCTATGCTATCCGATGCTCGAGTCGCTGAGTATGC	5'	Reverse Compliment



Polymerase Chain Reaction Quiz!



Polymerase Chain Reaction

5'	ACTGACGATACGATAGGCTACGAGCTCAGCGACTCATACG	3'	DNA Sequence
	3' ATGCTATC 5'		Reverse Primer
		5' AGCGACTC 3'	Forward Primer
3'	TGACTGCTATGCTATCCGATGCTCGAGTCGCTGAGTATGC	5'	Reverse Compliment



What product would we see?

Polymerase Chain Reaction

5'	ACTGACGATACGATAGGCTACGAGCTCAGCGACTCATACG	3'	DNA Sequence
	3' ATGCTATC 5'		Reverse Primer
		5' AGCGACTC 3'	Forward Primer
3'	TGACTGCTATGCTATCCGATGCTCGAGTCGCTGAGTATGC	5'	Reverse Compliment



What product would we see?


Polymerase Chain Reaction

5'	ACTGACGATACGATAGGCTACGAGCTCAGCGACTCATACG	3'	DNA Sequence
	3' ATGCTATC 5'		Reverse Primer
		5' AGCGACTC 3'	Forward Primer
3'	TGACTGCTATGCTATCCGATGCTCGAGTCGCTGAGTATGC	5'	Reverse Compliment



What product would we see?
NO PRODUCT. Why?

Polymerase Chain Reaction



What product would we see?
NO PRODUCT. Why?
Polymerase extension (5'→3')
would not yield geometric
amplification.

Polymerase Chain Reaction Quiz!

5'	ACTGACGATACGATAGGCTACGAGCTCAGCGACTCATACG	3'	DNA Sequence
		3' GCTGAGTA 5'	Reverse Primer
5'	GACGTTAC	3'	Forward Primer
3'	TGACTGCTATGCTATCCGATGCTCGAGTCGCTGAGTATGC	5'	Reverse Compliment




Polymerase Chain Reaction Quiz!

5'	ACTGACGATACGATAGGCTACGAGCTCAGCGACTCATACG	3'	DNA Sequence
		3' GCTGAGTA 5'	Reverse Primer
5'	GACGTTAC	3'	Forward Primer
3'	TGACTGCTATGCTATCCGATGCTCGAGTCGCTGAGTATGC	5'	Reverse Compliment



What product would we see?
Same product as before mostly.
What potential issue is here?

Polymerase Chain Reaction Quiz!



What product would we see?
Same product as before mostly.
What potential issue is here?
Physically, mismatch causes bubble,
but probably no issues.

Polymerase Chain Reaction Quiz!

5' ACTGACGATACGATAGGCTACGAGCTCAGCGACTCATACG 3' DNA Sequence
3' CTGAGTA 5' Reverse Primer
5' GACGATAC 3' Forward Primer
3' TGACTGCTATGCTATCCGATGCTCGAGTCGCTGAGTATGC 5' Reverse Complement



What product would we see?

Polymerase Chain Reaction Quiz!

5' ACTGACGATACGATAGGCTACGAGCTCAGCGACTCATACG 3' DNA Sequence
3' CTGAGTA 5' Reverse Primer
5' GACGATAC 3' Forward Primer
3' TGACTGCTATGCTATCCGATGCTCGAGTCGCTGAGTATGC 5' Reverse Complement



What product would we see?
No product. Polymerase can't bind
correctly to induce extension.

Polymerase Chain Reaction Quiz!

5' ACTGACGATACGATAGGCTACGAGCTCAGCGACTCATACG 3' DNA Sequence
3' CTGAGTA 5' Reverse Primer
5' GACGATAC 3' Forward Primer
3' TGACTGCTATGCTATCCGATGCTCGAGTCGCTGAGTATGC 5' Reverse Complement

PCR

What product would we see?
No product. Polymerase can't bind correctly to induce extension.

For the purposes of this assignment, ignore issues with bumps and first base mismatch.

PCR for Bacteria Identification

Do I want to identify all species of bacteria present in a sample (next gen sequencing) or determine whether a single species is present?

Both require PCR...


Identification of all bacteria present...

For sequencing, we want to make lots of copies of DNA for all bacteria present. Then we can sequence these.

The best way to do this is to look for a common gene across all bacteria and use PCR to make copies of that gene.

16s ribosomal RNA subunit


- Portion of ribosome (purple).
- Ribosomal function is essential for life so the sequence is mostly conserved across bacterial species.
- But we still see mutations in this gene. Why?



Task 3

How would we design our PCR primers for sequencing?

- We want to make copies of DNA in the 16s rRNA gene regardless of the species.



How would we design our PCR primers for sequencing?


- We want to make copies of DNA in the 16s rRNA gene regardless of the species.



What parts of our gene allow us to identify species?

How would we design our PCR primers for sequencing?

- We want to make copies of DNA in the 16s rRNA gene regardless of the species.




What parts of our gene allow us to identify species?

Where should we design our primers to bind such that when we sequence the PCR product, we can identify any species?

- Generate a function which takes as input a list of sequences (strings) and returns a forward and reverse primer which will generate a PCR product for sequences in the list.

Task 4

How would we design our PCR primers to detect presence of each of three bacterial species?



- Generate a function which takes as input a list of sequences (strings) and returns a set of primer pairs which will generate products for each sequence specifically (or *None* if the task is impossible).

Strategy Hints:

n = number of sequences in list

- 2 primers:
- $n + 1$ primers:
- $2n$ primers:

Task 5

- For your group's bacterial sequences, run Task 3 and Task 4. Generate a text file where each line contains an identifier and the sequence separated by a space.

- Example:

Group_10_JK_Task1_FWD ACTGCTACGGACGACT

Group_10_JK_Task1_REV TCAGCGACGAACGCTCT

Task 6

- Describe the design of an experiment to determine whether or not any of your n bacteria in Task 5 are present in a random sample of DNA. Be sure to include the list of reactions including the template DNA and primers used for each.